

PAT Basic® - Affordable Pediatric Auscultation Manikin, light skin tone

PRODUCT ID 1020098
WEIGHT 4.5 kg
BRAND Cardionics USA



KEY FEATURES

- The SimScope stethoscope is used to simulate pediatric heart, lung, bruit, and bowel sounds, utilizing a large sound library incorporated into the SimScope stethoscope.
 - The SimScope Wi-Fi enables wireless communication between the SimScope and the PAT Basic auscultation manikin, allowing users to change sounds remotely or to use in conjunction with a scenario-based simulation.
 - By placing the SimScope over the anatomically correct locations of PAT Basic, the programmed sounds will be played through the scope, allowing students to employ and enhance their diagnostic skills.
 - PAT Basic can be used for:
 - Individual learning, Group lectures, Viewing phonocardiograms.
 - Cardionics auscultation trainers now include a COVID-19 lung sound!
 - How is PAT Basic Student Auscultation Manikin different from other pediatric auscultation products available on the market?
 - It is lightweight & portable
 - It utilizes the extensive Cardionics proprietary pediatric heart, lung, bruit and bowel sounds library
 - Students can listen at anatomically correct auscultation sites
 - A customized library is available to fit specific scenarios and programs (contact us to find out more.)
- The SimScope Wi-Fi allows for remote selection of sound scenarios

PRODUCT DETAILS

The PAT BASIC trainer engages the technology of the SimScope stethoscope to simulate pediatric heart, lung, and bowel sounds utilizing a large sounds library. Using the provided SimScope stethoscope, users are able to access a large variety of pathological conditions in order to customize scenarios.

By placing the SimScope over the anatomically correct locations of PAT BASIC, the programmed sounds will be played through the scope, allowing for students to employ and enhance their diagnostic skills.

The optional SimScope WiFi upgrade enables wireless communication between the Simscope and PAT BASIC manikin, allowing users to change sounds remotely or to use in conjunction with a scenario-based simulation. PAT BASIC consists of a male manikin, a SimScope or SimScope WiFi, the SimScope software and site license, a USB cable, 1 SAM t-shirt and a user's manual.

