

## SAM Basic® - Affordable Adult Auscultation Manikin, light skin tone

PRODUCT ID 1020097  
WEIGHT 6.6 kg  
BRAND Cardionics USA



### KEY FEATURES

- The SimScope stethoscope is used to simulate adult heart, lung and bowel sounds, utilizing a large sound library incorporated into the SimScope stethoscope.
- By placing the SimScope over the anatomically correct locations of SAM Basic, the programmed sounds will be played through the training stethoscope, allowing students to employ and enhance their diagnostic skills. SAM Basic can be used for:
  - Individual learning, Group lectures, Viewing phonocardiograms.
  - Cardionics auscultation trainers now include a COVID-19 lung sound!
  - How is SAM Basic Student Auscultation Manikin different from other auscultation products available on the market?
  - It is lightweight & portable
  - It utilizes the extensive Cardionics proprietary heart, lung, bowel and bruit sounds library
  - Students can listen at anatomically correct auscultation sites
  - A customized library is available to fit specific scenarios and programs (contact us to find out more)
  - The SimScope Wi-Fi allows for remote selection of sound scenarios.
  - The SimScope Wi-Fi enables wireless communication between the SimScope and the SAM Basic manikin, allowing users to change sounds remotely or to use in conjunction with a scenario-based simulation.

### PRODUCT DETAILS

Auscultation is an essential clinical skill needed to assess and monitor patients' conditions. With Cardionics' SAM Basic Student Auscultation Manikin, teaching and learning this skill becomes even more versatile and economic. SAM Basic trainer can be used with the SimScope stethoscope (included). All sounds and videos are recorded from adult live patients for students to experience a life-like simulation.

SAM Basic Student Auscultation Manikin used in teaching and learning heart, lung and bowel sounds. The

---

SimScope stethoscope is used to simulate adult heart, lung and bowel sounds, utilizing a large sound library incorporated into the SimScope stethoscope. The SimScope WiFi upgrade enables wireless communication between the SimScope and the SAM Basic manikin, allowing users to change sounds remotely or to use in conjunction with a scenario-based simulation

By placing the SimScope over the anatomically correct locations of SAM Basic, the programmed sounds will be played through the scope, allowing for students to employ and enhance their diagnostic skills. SAM Basic can be used for:

- Individual learning
- Group lectures
- Viewing phonocardiograms

Delivery content:

- A portable SAM Basic Male Auscultation Torso
- One SimScope stethoscope
- The SimScope software and site license
- USB cable
- One SAM t-shirt
- User's manual

How is SAM Basic Student Auscultation Manikin different from other auscultation products available on the market?

- It is lightweight & portable
- It utilizes the extensive Cardionics proprietary heart, lung, bowel and bruit sounds library
- Students can listen at anatomically correct auscultation sites
- A Customized Library is available to fit specific scenarios and programs (contact us to find out more!) The SimScope WiFi upgrade (available separately) allows for remote selection of sound scenarios

Request an online demo today and find out how your students can benefit from having SAM Basic Trainer at hand!

All SAM models are available as light skinned or dark skinned version.